

# REBOOTS THERMO SLEEVE

HEAT AND COLD  
THERAPY

Thermo Sleeves can be used for heat and cold therapy, providing targeted treatment for sore muscles and support an optimal recovery. The high-quality, reusable gel pad in soft polyester provides 360-degree heat or cold treatment while flexibly adapting to the respective body part. The three available sizes ensure a snug fit, allowing the sleeves to apply a slight compression to the treated area. Both the heat and cold treatments also combine well with the compression massage of the Reboots Recovery Boots and Pants to further support the effect of the compression massage.



## FEATURES

S, M, L

Available sizes



360-degree Heat or Cooling



Compatible with Reboots



Usable for legs and arms



Reusable heat/cold pack in recyclable zipper bag

### HIGH WEARING COMFORT

The Thermo Sleeves are made of a high quality gel pad in soft polyester. They flexibly adapt to the part of the body to which the heat or cold therapy is applied, exerting a slight compression. At the same time, they are non-slip and easy to clean.



### HEAT & COLD

The Thermo Sleeves can be heated or cooled very easily and can therefore be used for many applications. The heating is done via the microwave or the oven, the cooling via the refrigerator.



### THE PERFECT ADD-ON FOR YOUR REBOOTS

Both heat and cold treatment can be perfectly combined with the compression massage of the Reboots Recovery Boots and Pants to further support the compression massage.



## SIZE CHART

	S	M	L
Product dimensions	26cm x 33cm circumference	26cm x 33cm circumference	26cm x 53cm circumference
Suitable Curcumference	26 - 35cm	33 - 44 cm	50 - 63 cm

## APPLICATION AREAS

- Wrist & ankle
- Elbow
- Upper and lower arm
- Knee
- Thigh
- Calf



www.maht.gr

info@maht.gr

210 69 80 000

## THERMO SLEEVE FULL LEG

The Thermo Sleeve Full Leg is specially designed for use on the entire leg and seamlessly covers the area from the groin to the ankle. Thanks to its stretching ability, it fits snugly around the leg without slipping or restricting your freedom of movement. Like our other Thermo Sleeves, the Thermo Sleeve Full Leg can be used to apply cold and heat, providing targeted treatment for sore muscles, joints and injuries. The high-quality, reusable gel pad offers 360-degree heat or cold treatment and adapts flexibly to your leg. Both the heat and cold treatment can also be perfectly combined with the compression massage of the Reboots Recovery Boots and Pants to further support the effect of the compression massage.



### FEATURES



Covers the entire leg from groin to ankle



360-degree Heat or Cooling



Compatible with Reboots



Usable for the legs



Reusable heat/cold pack in recyclable zipper bag

### APPLICATION OVER THE ENTIRE LEG

The Thermo Sleeve Full Leg extends from the groin to the ankle, guaranteeing seamless cold or heat application over the entire leg. Encased in soft polyester, the high-quality gel pad ensures that the target temperature is transmitted 360 degrees to all areas of your leg.



### HEAT & COLD

The Thermo Sleeve Full Leg can be heated or cooled very easily and can therefore be used for many applications. The heating is done via the microwave or the oven, the cooling via the refrigerator.



### THE PERFECT ADD-ON FOR YOUR REBOOTS

Both heat and cold treatment can be perfectly combined with the compression massage of the Reboots Recovery Boots and Pants to further support the compression massage.



### SIZE

	One size
<b>Product dimensions</b>	Length: 75 cm; 34 cm circumference at the top, 16 cm circumference at the bottom
<b>Suitable Circumference</b>	30 - 44 cm at the top of the thigh

### APPLICATION AREAS

- Entire leg: Use is possible on both legs or on one leg only.
- Use in combination with Reboots